Instruction Manual

FLASH lite 1401™

LED Dental Curing Light

For Sales Orders and Product Information Contact:

Manufactured by:

Discus Dental (USA)
8550 Higuera Street
Culver City, CA 90232
Tel: 800 422 9448
info@discusdental.com
www.discusdental.com

Distributed by:

Canada
Tel: 800 278 8282
www.discusdental.com

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands
Tel: +31 654 266 449
Europe@emergogroup.com

Asia Pacific

Australia
Tel: +61 2 8845 6900
Toll Free: 1800 621 448
www.discusdental.com/au

New Zealand
Tel: +61 2 8845 6900
Toll Free: 0800 600 041
www.discusdental.com/nz

All Other International
Tel: 310 845 8260
international@discusdental.com

Or contact your local distributor for further information.

DISCUS™

19-2497 051111

Instruction Manual

FLASH lite 1401™

LED Dental Curing Light

Powerful
Cordless
Compact

Instruction Manual
Manual de Instrucciones
Manual de instruções
Manuel d'utilisation
Manuale d'istruzioni
Bedienungsanleitung
Gebruiksaanwijzing
Brugsanvisning
Bruksanvisning
Instrukcja obsługi
Kullanım Kılavuzu
Руководство по эксплуатации
使用手册
使用説明書
사용 설명서

FLASH^{light} 1401™

Instruction Manual

Contents / Contenido / Conteúdo / Contenu / Materiale in dotazione / Inhalt / Inhoud / Indhold / Innehåll / Zawartość / İçindekiler / Содержание / 目次 / 目录 / 내용

- English - FLASHlite 1401 Instructions Manual
Corresponding Illustrations
- Español - FLASHlite 1401 Manual de Instrucciones
Ilustraciones correspondientes
- Português - Manual de instruções do FLASHlite 1401
Ilustrações correspondentes
- Français - Manuel d'instructions de la lampe Flashlite 1401
Avec illustrations
- Italiano - Manuale di istruzioni Flashlite 1401
Illustrazioni corrispondenti
- Deutsch - FLASHlite 1401 Bedienungsanleitung
Dazu gehörende Abbildungen
- Nederlands- FLASHlite 1401 Gebruiksaanwijzing
Bijpassende illustraties
- Dansk - FLASHlite 1401 brugsanvisning
Tilhørende illustrationer
- Svenska - FLASHlite 1401 bruksanvisning
Motsvarande bilder
- Polski – Instrukcja obsługi FLASHlite 1401
Odpowiednie ilustracje
- Türkçe - FLASHlite 1401 Kullanım Kılavuzu Manual
İlgili Çizimler
- Русский - FLASHlite 1401 Руководство по эксплуатации
Соответствующие иллюстрации
- 中文 - FLASHlite 1401 使用手册
相应插图
- 日本語 - FLASHlite 1401 取扱説明書
対応する説明書
- 한국어 - FLASHlite 1401 사용설명서
당 그림 설명

Product Description

The Discus Dental FLASHlite 1401 (see Figure 1) is a revolutionary LED dental curing light utilizing a light emitting diode (LED) for the polymerization of composite materials that contain the photo-initiator camphorquinone (CQ). The incorporation of an advanced LED, micro-controller and a lithium-ion battery enables the FLASHlite 1401 to be one of the most advanced, portable and powerful LED dental curing lights available today, offering curing power equal to or better than standard halogen curing lights.*

The FLASHlite 1401 uses the most sophisticated LED technology to provide the required wavelength range (460-480 nanometers) needed to ensure quick, reliable and efficient polymerization of products made with camphorquinone.

The FLASHlite 1401 is designed to be incredibly powerful yet easy to use. There is no need to memorize elaborate button sequences in order to determine cure times. The FLASHlite 1401 uses a one-touch operation for on/off and features an audible indicator for curing times at 10 second intervals. The light-weight, ergonomic design, featuring a curved extension, allows for maximum exposure to all areas of the oral cavity for excellent reach and super-efficient curing.

* Study available upon request.

Contents

(See Fig. 1)

- AC/DC Wall Plug Adapter
- Charger Base with Built-in Radiometer
- FLASHlite 1401
- One Replacement Lens Cap
- One Triangular Multi-tiered Hardness Curing Disk
- Instruction Manual

Technical Data

Battery:	Lithium-Ion 3.6V nominal 1800mAh
Wavelength Range:	460-480 nm
Light Intensity:	≥1100 mW/cm ²
Total Continuous Runtime with Fully Charged Battery:	25 minutes
Dimensions:	Diameter: 22 mm (.86 in) Length: 198 mm (7.8 in) 3.8oz
Weight:	3.8oz
AC/DC Wall Plug Adapter – Class II	
Input Voltage:	100-240V ~/ 50-60Hz / 400mA
Output Voltage:	9V=1.5A
Dimensions:	Height: 64.5 mm (2.5 in) Diameter: 62.7 mm (2.46 in)
Weight:	3.7 oz.
Part Number:	CR 1054

Charger Base with Built-in Radiometer – Class II
Time to Charge Discharged Battery: 4.5 hours
Operating Temperature: 15°C - 45°C (55°F - 104°F)
Total Height with Handpiece Inserted: 223.5 mm (8.8 in)
Handpiece – Internally Powered

Safety

1. The FLASHlite 1401 must be used in strict accordance with the following operating instructions. Discus Dental accepts no liability for any damage resulting from the use of this unit for any other purpose.
2. Prior to start-up of FLASHlite 1401 make sure that the operating voltage stated on the rating plate is compatible with the available mains voltage. Operation of the unit at a different voltage may damage the unit. Do not touch pins in charger base.
3. Exposure must be restricted to the area of the oral cavity in which clinical treatment is intended.
4. Operator must take care to avoid patient contact of the device while in use. Device generates heat while light is on. Device is not intended for direct patient contact. If contact is made, discomfort may occur.
5. Use only chargers and adapters supplied by Discus Dental. The use of third party chargers and/or AC/DC wall plug adapters will result in damage of the unit.
6. Do not allow cleaning agents to directly enter the unit as this may cause premature failure. **See Maintenance and Care.**
7. Only a Discus Dental authorized representative can service the unit. Please call Discus Dental for any questions. Phone numbers are on the back of this manual.
8. Do not use on patients with a history of photo biological conditions including urticarial solaris or erythropoetic protoporphyria or who are on photosensitizing medications.
9. Suitable blue-light filtering safety goggles must be worn by patient and dentist during use.

Classification

- Conforms to IEC 60601-1, EN60601-1, IEC 60601-1-2, CAN/CSA C22.2 No.601.1-M90

Unpacking / Setup

1. Insure all parts contained within package match the contents list.
2. Remove all pieces from the box.
3. Plug the AC/DC wall plug adapter into the mating connector on the back of the charger base.
4. Plug the AC/DC wall plug adapter connector into your most convenient AC outlet.
5. Place the FLASHlite 1401 into the charger base. The amber light on the charger base will indicate the FLASHlite 1401 is now charging. To insure highest level of performance, allow the FLASHlite 1401 to fully charge prior to first use. On average, it will take approximately 4-5 hours to initially charge FLASHlite 1401. Indicator will turn green when unit is fully charged. Keep in charger while not in use.

Operation

FLASHlite 1401 is supplied with the lens cap attached. A simple one-touch on/off button operation has been designed. Simply press the button once to activate and de-activate the light.

Duty Cycle

Intermittent use of 10 to 40 second curing times, with twice the cooling time, within maximum of 6 minutes. Place device in charger base to cool when not in use.

FLASHlite 1401 is equipped with a microprocessor to monitor the functionality of the unit including remaining battery life. FLASHlite 1401 is equipped with several audible signals including:

Timing Beeps

One beep at 10 seconds, two beeps at 20 seconds, three beeps at 30 seconds, four beeps at 40 seconds, five beeps

at 50 seconds, 6 beeps at 60 seconds and will immediately repeat this schedule continuously until manual shut-off by user by pressing the power button (no beeps).

Alarms

15 slow beeps—low battery warning

15 fast beeps—overheat alarm.

Note: When the external temperature reaches 50°C, the unit shuts off and then beeps 15 fast beeps.

Lens Cap (See Fig. 2)

The lens cap should be inspected prior to each use for blemishes, scratches, cracks or foreign substances that may impair the optical output of FLASHlite 1401.

It is recommended that lens caps be replaced every 2-4 weeks (Item #CR1024). To replace, simply unscrew the existing lens cap counter clockwise and replace by turning clockwise. Do not over tighten. An O-ring is mounted under the lens cap for alignment. Keep O-ring in place during lens cap replacement.

Barrier Sleeve

To provide optimum protection, we recommend you apply a clean plastic barrier sleeve for each new patient use. (Item # CR1033)

Cure Times

Due to the variation in VLC (visible light cured) materials, curing times will differ. Review manufacturer product instructions for recommended curing times. If using Matrixx by Discus Dental, these composite materials perform optimally when cured for approximately 20-40 seconds in 2-3 mm increments. **ALWAYS bench test new materials before use in vivo.**

Multi-Tiered Triangular Hardness Disk

The new multi-tiered triangular hardness disk enables you to measure various depths of composite shades. Each point of the triangle offers a different depth (2, 3 or 4 mm). The triangular disk is made of a plastic material, which mimics the properties of cured composite. To measure the

depth of cure, place the desired amount of composite into the appropriate tip well of the “cure” side of the triangle. Cure the material for the length of time suggested by the manufacturer. Next, check the hardness of the cured material by scraping the “test” side surface with a tungsten carbide instrument. Measure its hardness against that of the surrounding surface of the disk. The surface of the cured material should feel equal to or harder than that of the disk material. The disk has a Barcol hardness of 75 +/- 5 and represents a complete cure.

Note: If the cured composite material is softer than the surrounding disk surface, this may be the result of insufficient light output or compromised composite material. In order to verify light output, make sure the unit is fully charged and test intensity on the built-in radiometer on the FLASHlite 1401 charger base. Also review the manufacturer’s instructions for composite material.

Charger with Built in Radiometer

The FLASHlite 1401 charger is designed to maintain a peak level of curing performance. To charge, place the FLASHlite 1401 into the charger base contact side down. (See Fig. 3) The light indicator (See Fig 4) on the charger base will provide you with the information needed:

- Amber Light – Unit is charging
- Green Light – Unit is 100% Charged

Note: After initially fully charging the battery, FLASHlite 1401 can be used at any time during the charging process, even if the light is indicating a charging status (amber light).

For optimal performance, place FLASHlite 1401 back in charger base when not in use.

Measuring Output Intensity

The light intensity ($\geq 1100\text{mW/cm}^2$) can be measured for accuracy using the light indicator built into the charging base (See Fig. 4). To measure the intensity of the light being emitted from the unit, place the lens tip directly

onto the black aperture located on the charger base. Holding the FLASHlite 1401 directly onto the aperture, activate the FLASHlite 1401 by pressing the blue on/off button. Built within the charger base are 4 blue indicator lights. The indicator lights are not visible until the FLASHlite 1401 tip has been placed on the aperture and the FLASHlite 1401 turned on. The number of indicator lights that illuminate will identify the output of the light:

- **4 lights: 100%** of the light intensity available
- **3 lights: 80%** of the light intensity available
- **2 lights: 60%** of the light intensity available
- **1 light: 40%** of the light intensity available

The light intensity can be accurately determined only with the charger of the FLASHlite 1401. Do not use portable and/or handheld radiometers to measure the light intensity output of FLASHlite 1401 as these types of radiometers have been shown to have variability in accuracy and precision in measuring total spectral output. To measure output intensity of your FLASHlite 1401, only use the provided charger.

Always bench test new materials before use *in vivo*.

Maintenance and Care

Put FLASHlite back in charger after every use. Take care to keep charging contacts on the bottom of the device and charger pins in the charger base clean and dry. Do not use undo force against lens cap as this could damage threads and prevent proper replacement to lens cap. Do not take off the O-ring.

Recommended disinfecting agents, damp cloth:

- ≤ 3.4% Gluderaldehyde solution
- 70% Isopropyl alcohol
- 1.5% Phenol
- Quaternary ammonium

Take care to avoid wetting internal parts.

Disposal: Dispose of device with battery according to the corresponding legal requirements in your country. Please contact your local disposal contractor.

Transport and Storage


Ambient Temperature	4F to 104F (-20C to 40C)
Relative Humidity	10% to 90%
Atmospheric Pressure	.5atm to 1.0atm (500hPa to 1060hPa)

Warranty

Discus Dental warranties your FLASHlite 1401 against defects in material and workmanship for 24 months from date of purchase with proper usage. During that 24 month warranty period, Discus Dental will repair or replace a defective unit.

Please call Discus Dental Customer Service for any additional questions, comments or product information. Phone numbers are on the back of this manual.

Defects caused by misuse, neglect, accident, or abuse are not covered by this warranty.

Trouble Shooting Guide

	Problem	Causes	Solution
1	Continuous slow beeps	Low battery.	Recharge battery. Light will turn off when battery is too low to keep light on.
2	Continuous fast beeps	Device too warm.	Place hand held device in charger. Allow device to cool.
3	Three (3) Beeps, No light No light, no beeps	Device malfunction.	Contact customer service.
4	Light flickers / flashes	Low battery.	1. Recharge device for 3-5 days continuously. 2. If problem persists, contact customer service.
5	Lens cap cracked	Damage due to dropping. Damage due to cleaning material over-exposure. Over-tightening of lens cap. Wear and tear.	1. Carefully remove lens cap by unscrewing counter-clockwise. 2. Replace with new lens cap by tightening clockwise.
6	Light is completely charged (green indicator), but beeps while in use	Handle bottom or charger pins contaminated with debris. Low battery.	1. Unplug the charger base. 2. Wipe pins in the base to remove debris. Do not allow any liquid to drip into the base. 3. Check bottom of handle for debris. Wipe with dry cloth to remove debris. 4. Plug in charger base and place hand held device into base.

Electromagnetic Emissions

FL1401 is intended for use in the electromagnetic environment specified below. The customer or the user of the FL1401 should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The FL1401 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The FL1401 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

**Electromagnetic Immunity**

FL1401 is intended for use in the electromagnetic environment specified below. The customer or the user of the FL1401 should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U_T (>95% dip in U_T) 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	< 5% U_T (>95% dip in U_T) for 0, 5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the FL1401 requires continued operation during power mains interruptions, it is recommended that the FL1401 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the FL1401, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measure field strength in the location in which the FL1401 is used exceeds the applicable RF compliance level above, the FL1401 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the FL1401.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3.0 V/m.

NOTE 3: U_T is the a.c. mains voltage prior to application of the test level.

Separation distance between portable and mobile RF communication equipment and the FL1401

The FL1401 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the FL1401 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the FL1401 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0.1	0.38 m	0.38 m	0.73 m
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum putput power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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EUROPEAN DECLARATION OF CONFORMITY

We, Discus Dental, Inc. declare under our sole responsibility that the product Flashlite 1401, CR1050 to which this declaration relates is in conformity with the following Directive and standards:

- Medical Device Directive of 14 June 1993 (MDD 93/42/EEC)
- EN 60601-1:1990 + A1:1993 + A2:1995 + A13:1996
- EN 60601-1-2:2001
- EN 60825-1:2001: 1994 + A1: 2002 + A2: 2001

The Technical Construction File is maintained at Discus Dental, Inc.
8550 Higuera Street, Culver City, CA 90232, USA.

The authorized representative located within the Community is:

EMDAR B.V.
Ijsselburcht 3
P.o.b. 5486
6802 EL Arnhem
The Netherlands
Telephone: 0031 70 345 8570
Fax: 0031 70 346 7299

Date of issue: 11-13-06.

Place of issue: Culver City, CA, USA

Steven Ziemba, Vice President of Regulatory Affairs and Quality Assurance