Impression Material & Bite Registration Directions for Use

-
- SPLASH!® • MegaBite® • Vanilla Bite**
• Peppermint Snap Clear Bite • Milk Chocolate Bite

Directions for Use	2
Mode d'emploi.....	10

Product Description:

Splash!® Extra Lite Body (Orange) and Lite Body (Green)

ISO 4823; Type 3: Low Consistency

Splash!® Medium Body (Blue) and Heavy Body (Purple)

ISO 4823; Type 2: Medium Consistency

Splash!® Putty (Purple)

ISO 4823; Type 0: Putty Consistency

Vanilla Bite Registration (White)

No specification existing.

Peppermint Snap Clear Bite Registration (Clear)

No specification existing.

Milk Chocolate Bite Registration (Brown)

No specification existing.

MegaBite® Registration (Teal)

No specification existing.

Range of Indications:

Splash! Impression material and Discus Dental Bite Registration systems are addition-reaction silicone formulas. Splash! is ideal for dual viscosity impression techniques both one-step and two-step as well as for a variety of dental procedures including, but not limited to, crowns, bridges, inlays, partial and complete dentures, denture repairs, implants and for forming functional peripheries. Bite Registration materials are ideal for capturing accurate interocclusal and jaw relation records. With the exception of MegaBite, all our Bite Registration materials, can be used in a tray or as a trayless quadrant material. MegaBite should only be used without a tray.

Working & Setting Times

	Max. Working Time	Total Set Time
Splash! Viscosities	2 min. 15 sec.	4 min. 30 sec.
Splash! Putty & Putty Paks®	1 min. 15 sec.	5 min. 30 sec.
Splash! Half-Time Cartridge Viscosities	1 min. 0 sec.	2 min. 15 sec.
Splash! Half-Time Putty & Putty Paks®	1 min. 0 sec.	2 min. 45 sec.
Vanilla Bite Registration	20 sec.	55 sec.
Peppermint Snap Clear Bite Registration	30 sec.	2 min. 30 sec.
Milk Chocolate Bite Registration	45 sec.	1 min. 30 sec.
MegaBite Registration	20 sec.	55 sec.

Directions for Use

Instructions for Cartridge Viscosities:

1. Load the cartridge into the cartridge dispensing gun.
2. Remove cartridge closure and discard.
3. Before applying the mixing tip, extrude impression material until it exits both chambers at the same rate.
4. Attach the mixing tip on the cartridge and squeeze the dispenser to extrude the required amount of material directly onto the prepared tooth or into an impression tray or delivery syringe.
5. Leave the filled mixing tip on the cartridge after use.
6. Before using the cartridge again, remove and discard the used mixing tip. Before re-applying another mixing tip, be sure to extrude material so the product is exiting both chambers at the same rate.
7. Attach a fresh mixing tip and proceed as usual.
8. The cartridge dispensing gun can be cold sterilized.

Instructions for Putty Jars:

1. Scoop even amounts of base and catalyst from jars.
2. Make sure to keep base and catalyst scoops separate from each other.
3. Knead putty thoroughly until the color is homogenous and streak free.
Note: Do not mix while wearing latex gloves: vinyl or nitrile gloves are recommended.
4. Load mixed putty into tray.

Instructions for Putty Paks®:

1. Peel off lid from Putty Pak.
2. Fold base and catalyst together.
3. Squeeze product together.
4. Remove product from pak and mix thoroughly until color is homogenous and streak free.
Note: Do not mix while wearing latex gloves: vinyl or nitrile gloves are recommended.
5. Load mixed putty into tray.

Instructions for FatPak™ Viscosities:

Please refer to the illustrations in the back of this booklet.

1. Before first use, turn the white safety pin on the activating head of the large foil bag in the direction indicated by the arrow in the drawing, then remove and dispose of the pin (fig. 1).
2. Push both foil bags into the FatPak cartridge body. Be sure that the cartridge body and activating head notches are aligned. This will ensure optimum alignment of the outlet nozzles (fig. 2).
3. Press the activating heads onto the cartridge body until they lock firmly into final position (fig. 3 & 4). The foil bags are automatically pierced by a pin when the activating heads are pressed onto the cartridge body.
4. Insert the filled cartridge body into the automix unit.
5. Switch the device on and bleed about 1/4 inch of material to ensure even flow of base and catalyst.
6. Assemble the dynamic mixer tip and close the locking clip for the mixer located on the cartridge body. Note: use only Discus Dental FatPak dynamic mixer tips for optimum performance.
7. Fill the tray with the desired amount of material by holding the tray in a tilted position and gently pressing against the dynamic mixer. Keep the tip immersed in the material to avoid air bubbles (fig. 5).

8. Do not remove the used mixer tip from the cartridge body as it acts as a seal.
9. Before using again, release the locking clip to remove the used mixer and check that the outlets in the activating heads are not blocked. Fit a new mixer and close the locking clip and continue with steps 7 – 8.
10. To remove the cartridge from the automix unit, hold the top of the locking clip with one hand and remove the complete cartridge (fig. 6).
11. Remove and dispose empty foil bags including activating heads from the cartridge body.

Splash! & Bite Registration Technical Specifications:

- Cartridge Mixing Technique (XL, L): Use a 1:1/1:2 dispensing cartridge gun, with a yellow static mixing tip (4.2-12D Diameter).
- Cartridge Mixing Technique (M, H): Use a 1:1/1:2 dispensing cartridge gun, with a blue static mixing tip (6.5-11D Diameter), or 5:1 automix unit and green dynamic mixing tip for FatPak™ foil bags.
- Storage Conditions: In a dry place at 15 – 25 °C or 59 – 78 °F and relative humidity of 50% or less. MegaBite should be stored at 18 – 25 °C or 64 – 78 °F. Do not expose to direct heat or sun. Refrigeration will not extend shelf life. Always store used cartridges with mixing tips intact to protect remaining material in cartridges from contamination. Re-using the original cap could cause cross-contamination of base and catalyst components and result in premature setting or hardening of the material at the tip. Store cartridges upright.
- Disinfection: Use a cold sterilant.
- Electroplating: Impressions may be copper or silver plated.
- Casting: The impressions can be poured as soon as 30 minutes after removal from the mouth. Surface tension will be reduced if the impression is briefly rinsed out with a surfactant such as Discus Dental Speed Spray.
- Shelf Life: Splash! & Splash! Half-Time have a shelf life of three years from the date of manufacture. Bite Registrations, with the exception of MegaBite, have a shelf life of three years from the date of manufacture. MegaBite has a shelf life of two years from the date of manufacture.

Splash! & Splash! Half-Time Extra Lite Body

Color: Orange

Recovery from Deformation: 99.7%

Strain in Comp.: 4.5%

Shore A Hardness: 45

Max. Dimensional Change: <0.1%

Splash! & Splash! Half-Time Lite Body

Color: Green

Recovery from Deformation: 99.7%

Strain in Comp.: 3.5%

Shore A Hardness: 55

Max. Dimensional Change: <0.1%

Splash! & Splash! Half-Time Medium Body Color: Blue

	Cartridge	FatPak
Recovery from Deformation:	99.5%	99.5%
Strain in Comp.:	2.5%	max 4%
Shore A Hardness:	64	56-62
Max. Dimensional Change:	<0.1%	<-0.25%

Splash! & Splash! Half-Time Heavy Body Color: Purple

	Cartridge	FatPak
Recovery from Deformation:	99.2%	99.3%
Strain in Comp.:	2.3%	<3.5%
Shore A Hardness:	67	57-63
Max. Dimensional Change:	<0.1%	<-0.25%

Splash! & Splash! Half-Time Putty Color: Purple

Recovery from Deformation: 99%

Strain in Comp.: 2.1%

Shore A Hardness: 72

Max. Dimensional Change: <0.1%

Vanilla Bite Registration Color: White

Recovery from Deformation: 99.7%

Strain in Comp.: 0.7%

Shore A Hardness: 90

Max. Dimensional Change: <0.1%

Peppermint Snap Clear Bite Registration Color: Clear

Recovery from Deformation: 99.6%

Strain in Comp.: 1.9%

Shore A Hardness: 63

Max. Dimensional Change: <0.1%

Milk Chocolate Bite Registration Color: Brown

Recovery from Deformation: 99.7%

Strain in Comp.: 0.7%

Shore A Hardness: 90

Max. Dimensional Change: <0.1%

MegaBite Registration Color: Teal

Recovery from Deformation: N/A

Strain in Comp.: N/A

Shore D Hardness: >40

Max. Dimensional Change: <-0.20%

Special Notes & Precautions:

- Cartridge material must be used with a 1:1/1:2 dispensing cartridge gun system. FatPak material can be used in a 5:1 compatible automix unit with dynamic mixer tips.
- The temperature of Splash!, Splash! Half-Time and Discus Dental Bite Registrations affects the working time. Ideally it should be kept at room temperature (15-25 °C or 59 – 78 °F). Even within those boundaries, changes in temperature will affect the working time: a 2°C (3.6 °F) temperature increase will shorten working time 10-15 seconds; conversely, a 2 °C (3.6 °F) temperature decrease will lengthen working time 10-15 seconds. You may refrigerate Splash!, Splash! Half-Time and Dental Bite Registration materials to increase working and setting times.
- Always “bleed” the cartridge before applying a fresh mixing tip. This will ensure even dispensing of fresh base and catalyst pastes.
- For even more precise application, intraoral tips are available.
- Solutions applied prior to impression taking, such as retraction solution, may impede the setting reaction of the impression. Therefore, material must be cleaned off thoroughly with water and dried.
- Do not leave remains of the impression material in the sulcus.
- The possibility of allergic reactions to the product in sensitive persons cannot be fully excluded. Possible reactions may include temporary numbness or tingling in the soft tissue.
- When taking an impression adjacent to an existing fixed bridge, make sure all undercuts are sufficiently blocked out.
- Some latex gloves or hand lotions could seriously inhibit the setting of vinyl polysiloxane impression materials. All intraoral areas touched with latex gloves should be rinsed and dried prior to impression taking.
- Avoid contaminating Splash!, Splash! Half-Time and Discus Dental Bite Registration impression materials with chemicals that could

inhibit the set (e.g. acrylic and methacrylate residues, latex rubber and sulphur compounds).

- Certain hemostatic agents may interfere with the set of Splash!, Splash! Half-Time and Discus Dental Bite Registration impression materials. Care must be taken to wash and dry contaminated areas before applying the syringe material. All intraoral areas touched by latex gloves should be rinsed and dried prior to impression taking.

The information provided by Discus Dental is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use.

For professional use only. U.S.A. Federal Law restricts this device to sale by or on the order of a dentist.

Contra-Indications

There are no known contra-indications provided when the product is used as directed. There are no known harmful reactions or secondary effects on patients and/or dental personnel using vinyl polysiloxane impression materials.

Compatibility

Splash!, Splash! Half-Time & Discus Dental Bite Registrations are addition-reaction type silicone impression materials and should not be combined with standard or condensation type silicone materials or polyether materials.

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