

MATERIAL SAFETY DATA SHEET

TRADE NAME:	HOLISTIC GERISTORE CONDITIONER
--------------------	---------------------------------------

CHEMICAL NAME: Aqueous HNO₃ Oxalate Solution

I. PHYSICAL DATA

BOILING POINT (°C):	NE
MELTING POINT:	N/A
SPECIFIC GRAVITY (H ₂ O=1):	1.0
VAPOR PRESSURE (20°C):	NE
VAPOR DENSITY (Air=1):	Less than 1
SOLUBILITY IN WATER:	Soluble
EVAPORATION RATE (Butyl Acetate=1):	Less than 1
APPEARANCE AND ODOR:	Tannish liquid

II. HAZARDOUS INGREDIENT

PRINCIPAL HAZARDOUS COMPOUND:	<u>TEL</u>	<u>VEL</u>
Aqueous HNO ₃ Oxalate Solution	NE	NE

III. FIRE AND EXPLOSION HAZARD DATA

FLASH POINT:	N/A
EXTINGUISHING MEDIA:	N/A
SPECIAL FIRE EXTINGUISHING PROCEDURES:	N/A
UNUSUAL FIRE AND EXPLOSION HAZARDS:	N/A
LEL EXPLOSIVE LIMIT:	NE
VEL EXPLOSIVE LIMIT:	NE

IV. HEALTH HAZARD DATA / FIRST AID PROCEDURES

INHALATION:	None known.
SKIN:	Flush with soap and water.
EYES:	Flush with water for 15 minutes and consult physician.
INGESTION:	Take glass of milk and consult physician.
CARCINOGENICITY:	None known.

V. REACTIVITY DATA

STABILITY:	Stable
INCOMPATIBILITY:	Bases
HAZARDOUS DECOMPOSITION:	May form oxides of nitrogen
HAZARDOUS POLYMERIZATION:	None

VI. SPILL OR LEAK PROCEDURES

STEPS TO BE TAKEN IF MATERIAL IS RELEASED OR SPILLED:	Wipe with absorbent towel.
WASTE DISPOSAL METHOD:	Dispose according to government regulations.

EMERGENCY PHONE NUMBER: (805) 922-8491

While Den-Mat Corporation believes that the data contained herein are factual and the opinions expressed are those of qualified experts regarding the results of the tests conducted, the data are not to be taken as a warranty or representation for which Den-Mat assumes legal responsibility. They are offered solely for your consideration, investigation and verification. Any use of these data and information must be determined by the user to be in accordance with applicable Federal, State and local laws and regulations.

DEN-MAT CORPORATION, 2727 SKYWAY DRIVE, SANTA MARIA, CA 93455