

Dentists are raving about this new system!

LUMINEERS® Placement System



APPLICATIONS:

- PRE-TREATMENT
- SHADE SELECTION
- COLOR MODIFICATION
- FINISHING



All the tools you need for LUMINEERS placement in convenient, reusable containers.

The LUMINEERS Placement System has everything you need to place CERINATE® LUMINEERS. This kit eliminates a cabinet of four to six different kits for your LUMINEERS placement. No more trying to remember everything you need and having it ready at chairside. In a convenient, ergonomically designed package, the LUMINEERS Placement System takes you through impressions, try-in, placement to finishing for that final LUMINEERS smile.

In many cases, without any tooth modification, you can make an accurate impression with LUMINEERS Impression material and send it off to the Cerinate Smile Design Studios to make your patient's LUMINEERS. With the Ultra-Bond® Try-in Paste, you can decide with your patient the final shade for that new LUMINEERS smile. This will save you time and money on second appointments. Once you have decided on the final shade, everything is at chairside, right at your fingertips, for LUMINEERS placement. Pretreatment of LUMINEERS with Porcelain Conditioner,

Cerinate Prime and Tenure® S is quick and easy. The teeth to be bonded can then be etched and bonded with clinically proven Etch 'N' Seal and Tenure MPB. You are now ready to place that new LUMINEERS smile with Ultra-Bond Plus resin cement. For those cases where you want to create a polychromatic result, or for darkened, stained teeth, there is a pallet of Cerinate Shade Modifiers and TetraPaques. When you need to mask discolorations, the clinician's choice is TetraPaque®. In one thin application it can mask even the darkest tooth discolorations and stains. Once the LUMINEERS are in place, finishing and polishing is easy using the CeriSaw, CeriSander and the LUMINEERS finishing diamonds and burs.

The LUMINEERS Placement System is well-designed and uncomplicated. By using the components of the LUMINEERS Placement System, you are assured of the LUMINEERS 5-year limited warranty. Why would you use anything else?

“There is no product on the market today that rivals the Cerinate LUMINEERS. Beautiful, easy to work with and excellent fit. It has changed the personalities and brightened the smiles of all patients that use them.”

—DR. HOBART FOD, NEWPORT, TN

NO SHOTS, NO POST-OP SENSITIVITY



BEFORE



AFTER

KIT CONTENTS (PART #031381100):

Tray 1

- 1 Bottle Cerinate Prime® 6 g
- 1 Bottle Cerinate® Conditioner 6 g
- 1 Bottle Dry Bond 10 mL
- 1 Syringe Etch 'N' Seal 3 mL
- 1 Syringe PorceLock® Porcelain Etchant 1.75 mL
- 1 Syringe Paint-On Dental Dam 3 mL
- 1 Syringe Cerinate Shade Modifier White 1 g
- 1 Syringe Cerinate Shade Modifier Yellow 1 g
- 1 Syringe Cerinate Shade Modifier Blue 1 g
- 1 Syringe Cerinate Shade Modifier Violet 1 g
- 1 Syringe Cerinate TetraPaque® Dentin 1 g
- 1 Jar Cerinate TetraPaque Indicator 4 g
- 1 Syringe Cerinate TetraPaque Enamel 1 g
- 1 Syringe Cerinate TetraPaque Pink 1 g
- 1 Syringe Cerinate TetraPaque Violet 1 g
- 1 Bag of Light Curing Test Rings
- 1 Brush
- 60 Tips for Clean-up Brush
- 4 Dappen Dishes
- 1 Mixing Pad
- 6 Pink Tips for Etch 'N' Seal and PorceLock Syringes
- 20 Black Tips for Modifier Syringes

Tray 2

- 1 Syringe Ultra-Bond® Try-In Paste A1 1 g
- 1 Syringe Ultra-Bond Try-In Paste B1 1 g
- 1 Syringe Ultra-Bond Try-In Paste White 1 g
- 1 Syringe Ultra-Bond Try-In Paste B0 1 g
- 1 Jar Tenure® S Dab-Eze® 3 g
- 1 Syringe Ultra-Bond Plus Resin Cement A1 10 g
- 1 Syringe Ultra-Bond Plus Resin Cement B1 10 g
- 1 Syringe Ultra-Bond Plus Resin Cement White 10 g
- 1 Syringe Ultra-Bond Plus Resin Cement B0 10 g
- 1 Bottle Tenure A Adhesive 6 mL
- 1 Bottle Tenure B Adhesive 6 mL
- 1 Mixing Pad
- 20 Brown Auto-Mixing Tips for Resin Cement Syringes
- 12 Pink Tips for Try-In Syringes

Tray 3

- 1 Lumineers® Finishing Diamond and Bur Block – (4 Diamonds; 1 Bur)
- 1 Schure 349 Instrument
- 1 Syringe Lumineers Impression PVS Fast Set/Heavy Body 50 mL
- 1 Syringe Lumineers Impression PVS Fast Set/Light Body 50 mL
- 1 Syringe Porcelain Polishing Paste 5 g
- 8 Auto-mixing Tips for PVS Syringes
- 10 Pink Tips for Polishing Syringes
- 2 CeriSaw and CeriSander Handles
- 10 Ceri-Saw Anterior Blades
- 6 CeriSander Fine Strips
- 6 CeriSander Ultra-fine Strips



LIVE
LUMINEERS®
PLACEMENT
VIA LUMICAST™*

SEE A LIVE LUMINEERS® PLACEMENT
All Cerinate® Smile Design Workshops now feature a live LUMINEERS placement via LUMICAST™* —the pain-free alternative to traditional veneers. Call 1-800-445-0345 today to find a course near you.



*Subject to technical capabilities.

LUMINEERS® Placement System

Do you have additional questions? No problem, our trained professional representatives are here to help.
Call 1-800-445-0345 to speak to someone live today! Or visit us online at www.denmat.com.

DenMat®
ISO 9001 CERTIFIED

©2005 Den-Mat Corporation. World Rights Reserved. 821219600 11/05PAC