

Safety Data Sheet

according to Regulation (EU) No. 1907/2006

Moxie Total Etch Adhesive, light-cure

Print date: 09.21.2010

Page 1 of 4

1. Identification of the substance/preparation and of the company/undertaking

MSDS-ID: Total Etch Adhesive, light-cure [REF2050]

Company/undertaking identification

Company name: Discus Dental, LLC
Street: 8550 Higuera St
City: Los Angeles, CA 90232
Telephone: (800) 422-9448 Telefax: (310) 845-1537
Internet: www.discusdental.com

2. Hazards identification

Classification

Indications of danger : Irritant
R-phrases:
Irritating to eyes, respiratory system and skin.
May cause sensitization by skin contact.

3. Composition/information on ingredients

Chemical characterization (Mixture)
chemical characterization (preparation): Acrylate.-resin.

4. First aid measures

After inhalation

Move victim to fresh air. Put victim at rest and keep warm.

After contact with skin

After contact with skin, wash immediately with: Water and soap.

After contact with eyes

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

After ingestion

Immediately get medical attention.

5. Fire-fighting measures

Suitable extinguishing media

Water fog. Extinguishing powder. Sand. Foam. Carbon dioxide (CO₂).

Extinguishing media which must not be used for safety reasons

High power water jet.

6. Accidental release measures

Personal precautions

Wear suitable protective clothing. Provide adequate ventilation.

Environmental precautions

Do not empty into drains or the aquatic environment.

Methods for cleaning up/taking up

Collect mechanically.
Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents).

Safety Data Sheet

according to Regulation (EU) No. 1907/2006

Moxie Total Etch Adhesive, light-cure

Print date: 09.21.2010

Page 2 of 4

7. Handling and storage

Handling

Advice on safe handling

Keep container tightly closed. Wear suitable protective clothing and gloves. Avoid contact with eyes.

Storage

Requirements for storage rooms and vessels

Keep container tightly closed in a cool, well-ventilated place.
Store only in original container.

8. Exposure controls/personal protection

Exposure limit values

Exposure controls



Occupational exposure controls

Ensure adequate ventilation of the storage area.

Protective and hygiene measures

When using do not eat or drink.

Respiratory protection

The following must be prevented: inhalation.

Hand protection

Tested protective gloves are to be worn: Suitable material: NBR (Nitrile rubber).

Eye protection

Tightly sealed safety glasses.

9. Physical and chemical properties

General information

Physical state: liquid
Colour: yellow
Odour: like: ester

Important health, safety and environmental information

	Test method
pH-Value (at 20 °C):	3,5
Changes in the physical state	
Boiling point:	100 °C
point of decomposition:	> 200 °C
Flash point:	> 100 °C
Vapour pressure: (at 20 °C)	23 hPa
Density:	1,1 g/cm ³
Solubility in other solvents:	Solubility ca. 50 %

Safety Data Sheet

according to Regulation (EU) No. 1907/2006

Moxie Total Etch Adhesive, light-cure

Print date: 09.21.2010

Page 3 of 4

Vapour density: > 1

10. Stability and reactivity

Conditions to avoid

Light. heat.

Materials to avoid

Keep away from strong acids, leachates, heavy metal salts and reducing materials.

Hazardous decomposition products

Can be released in case of fire: Gas / vapours, irritant. (Acrylate., pungent)

11. Toxicological information

Acute toxicity

LD50: Rat. 5500mg/kg

Corrosive and irritant effects

Frequent or prolonged contact with skin may cause dermal irritation.

Irritant effect on the eye: Frequent or prolonged contact with eye may cause irritation.

Sensitising effects

May cause sensitization by skin contact.

Additional information on tests

Contains Methacrylates: May produce an allergic reaction. May cause sensitization by skin contact.

12. Ecological information

Ecotoxicity

Aquatoxicity LC50: > 360 mg / l / 48h

Persistence and degradability

Product is partially biodegradable. 50%

Further information

Do not empty into drains or the aquatic environment. Leakage into the environment must be prevented.

13. Disposal considerations

Advice on disposal

Can be burnt together with household waste in compliance with official regulations in contact with approved waste disposal companies and with authorities in charge.

Carry out a burning of hazardous waste according to official regulations.

Waste disposal number of waste from residues/unused products

180106 WASTES FROM HUMAN OR ANIMAL HEALTH CARE AND/OR RELATED RESEARCH (except kitchen and restaurant wastes not arising from immediate health care); wastes from natal care, diagnosis, treatment or prevention of disease in humans; chemicals consisting of or containing dangerous substances
Classified as hazardous waste.

14. Transport information

Other applicable information

Not a hazardous material with respect to these transportation regulations.

15. Regulatory information

Safety Data Sheet

according to Regulation (EU) No. 1907/2006

Total Etch Adhesive, light-cure [REF2050]

Print date: 09.07.2010

Page 4 of 4

Labelling

Danger symbols: Xi - Irritant



Xi - Irritant

Hazardous components which must be listed on the label

Methacrylic esters.

R phrases

Irritating to eyes, respiratory system and skin.
May cause sensitization by skin contact.

S phrases

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
After contact with skin, wash immediately with plenty of Water and soap.

National regulatory information

Water contaminating class (D): 2 - water contaminating

16. Other information

Full text of R-phrases referred to under sections 2 and 3

Irritating to eyes, respiratory system and skin.
May cause sensitization by skin contact.

Further Information

-

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)