

ALL ABOUT THE CURVES



**SOL**<sup>TM</sup>  
Power of light

Whether you use a laser for gingival contouring, tissue retraction for impressions or laser-assisted periodontal therapy, the SOL laser is the ideal tool for soft-tissue management.



### Troughing

Tired of packing a retraction cord? Laser troughing creates separation between the hard tissue margins and the gingiva while simultaneously stopping any tissue bleeding prior to taking an impression or restoring a sub-gingival margin.

### Gingivectomy

Unwanted tissue can be removed within a bloodless field, with minimal post-operative pain. The result is an increased crown length for improved restoration retention and a symmetrical smile line. Final impressions for restorations can be completed in the same appointment, saving you time and money.



**SOL. Everything you want in a laser.**

Call DenMat at 800.445.0345 or visit [DenMat.com](http://DenMat.com).




# INTRODUCING SOL™

THE NEW LASER TECHNOLOGY DESIGNED  
FOR THE WAY YOU PRACTICE.

SOL™ delivers true portability, power and precision, along with enhanced ergonomics, simple interface, and a high-contrast aiming beam designed to make laser dentistry easier than ever. It also happens to be incredibly affordable.

Great curves. Outstanding performance. Unbeatable value.  
SOL is the most attractive desktop laser in dentistry today...  
and it's from DenMat.



Lasers have revolutionized the way today's dental professionals practice. Compared to traditional methods, diode lasers perform a wide array of procedures with far less collateral tissue damage and immediate hemostasis, promoting faster healing with less gingival recession. Your patients benefit from a higher standard of care, with faster healing times and greater comfort without the need for local anesthetic.

SOL is the ideal device for both entry-level and experienced laser users, delivering all the power and essential features you demand in a laser, but a fraction of the price. SOL combines true portability and superior precision at a value that no other laser can match. Whether you're upgrading your arsenal, arming your Hygiene team, or simply seeking your first diode, look no further than SOL.

### TRUE PORTABILITY


Internal battery allows you to move between operatories for over 3 hours without plugging in, and the contained fiber spool system keeps unwanted cord length out of your way.



Battery indicator

### SIMPLICITY


3 easy-to-use settings designed to simplify your treatment protocols, with a programmable Custom setting to address your unique needs.



Programmable presets

### VISIBILITY


Blue aiming beam enhances contrast on tissue so you know where's you're working.



Blue aiming beam

### ERGONOMICS

Curved handpiece design provides better access to the posterior in multiple positions and reduces fatigue during longer procedures.



Curved handpiece design

### SINGLE-USE

Disposable tips eliminate the hassle and inconsistency of strip-and-cleave systems.



Disposable tips